

Alpers Lab

In Situ Hybridization

Use RNase free glassware, containers, slide carriers and solutions. Use RNaseZap to rinse containers, etc and ddH₂O (or DEPC treated dH₂O). Wear gloves!

Solutions - Day 1 (RNase free)

0.5x SSC
Proteinase K
Proteinase K buffer
Prehybridization solution
Box buffer

Solutions - Day 2

2x SSC
RNase A
RNase buffer
0.1x SSC + 0.5% Tween 20

DAY ONE

Use only green labelled (RNase free) solutions, glassware, etc.

1. Warm 200 ml Proteinase K buffer in 37°C water bath. Thaw proteinase K aliquot.
2. Label all slides with specimen number, probe and number of days to be exposed (or #1 and #2)
3. Deparaffinize:

Xylene x 3	10 min. each
100% Ethanol x 3	1-2 min. each
95% Ethanol x 2	1-2 min. each
70% Ethanol x 1	1-2 min. each
0.5x SSC rinse x 3	10 secs. each
RNase buffer rinse	10 secs. each

****Use clean xylenes, alcohols****

4. Add 100 ul stock proteinase K into 200 ml Proteinase K buffer just prior to adding slides. Incubate slides at 37°C for 30 minutes. (Note: the concentration of proteinase K can be varied to optimize hybridization - depends on type of tissue being used, but 5 ug/ml works well for kidney)
5. Wash 0.5xSSC x 3 30 secs each

Optional - Dehydrate slides using graded ethanols and allow to air dry

6. Make sure oven is warmed to 50°C. Lightly saturate paper lining the incubation boxes with box buffer. Boil the prehybe DTT solution for 5 minutes.
7. Blot the tissue sections on filter paper (not necessary if dehydrated and air dried), add 100 ul of prehyb DTT to the tissue. Place in incubation boxes and cover tightly.
8. Incubate for 2 hours at 50°C.
9. Determine the amount of probe needed - plan for 50 ul of hybridization solution per slide (make extra as solution is very viscous). Add approximately 500,000 cpm per slide.

$$(50 \text{ ul}) * (\# \text{ of slides}) = \text{ul prehyb. solution needed}$$

$$500,000/(\text{cpm probe}) = \text{ul probe per slide}$$

$$(\text{ul probe per slide}) * (\# \text{ of slides}) = \text{ul probe to add to hyb. solution}$$

10. Transfer the volume of probe needed to a microcentrifuge tube and denature the probe by placing in a 95°C water bath for 3 minutes.
11. Place the probe in an ice bath.
12. Mix the required amount of prehybridization solution with the probe. Mix well. This is the **hybridization solution**.
13. After the slides have been prehybridized for 2 hours, add 50 ul of the hybridization solution directly to the slides.
14. Incubate overnight at 50°C in humidity boxes.

